

ChimeraLM detects amplification artifacts for accurate structural variant calling in long-read single-cell sequencing

Yangyang Li^{1†}, Qingxiang Guo^{1†}, Rendong Yang^{1,2*}

¹Department of Urology, Northwestern University Feinberg School of Medicine, 303 E Superior St, Chicago, 60611, IL, USA.

²Robert H. Lurie Comprehensive Cancer Center, Northwestern University Feinberg School of Medicine, 675 N St Clair St, Chicago, 60611, IL, USA.

*Corresponding author(s). E-mail(s): rendong.yang@northwestern.edu;

Contributing authors: yangyang.li@northwestern.edu;

qingxiang.guo@northwestern.edu;

[†]These authors contributed equally to this work.

Abstract

Single-cell genomics enables unprecedented insights into cellular heterogeneity across cancer, development, and neuroscience, but faces a fundamental challenge: Whole Genome Amplification (WGA) introduces chimeric artifacts that generate false Structural Variations (SVs), undermining biological interpretations. Current computational methods fail to distinguish these amplification-induced artifacts from genuine genomic rearrangements. Here we present ChimeraLM, a genomic language model that resolves this bottleneck by learning sequence-level features that discriminate authentic biological sequences from WGA artifacts. Built on Hyena operators for efficient long-range modeling, ChimeraLM analyzes nanopore sequencing reads at single-nucleotide resolution, achieving 95% recall with 70% precision on held-out test data. When applied to WGA-amplified cells, ChimeraLM reduces chimeric read content by ~90% (from 46% to 4.9% on PromethION, 23% to 1.5% on MinION) while preserving 87–92% of true SVs. This filtering improves SV validation rates 10–16 fold and eliminates the characteristic false-positive inversion (INV) bias inherent to unprocessed WGA data, restoring SV type distributions to bulk-like profiles. Attention weight visualization reveals that ChimeraLM focuses on chimeric junction regions with

single-base precision, learning biologically interpretable features of template-switching events. By enabling confident detection of chromosomal instability and structural variation at single-cell resolution, ChimeraLM addresses a fundamental data quality barrier across applications in cancer genomics, developmental biology, and precision medicine. The model is available at <https://github.com/ylab-hi/ChimeraLM>.

Keywords: Whole Genome Amplification, Single Cell, Genomic Language Model, Structural Variation

Main

Single-cell genomics has revolutionized our understanding of cellular heterogeneity by enabling characterization of individual cells rather than bulk populations [1–4], revealing previously hidden biological complexity. This approach has proven instrumental in uncovering rare cell types [4], tracking developmental trajectories [3], and elucidating tumor evolution through clonal architecture analysis. However, the limited DNA content in a single cell—typically only 6–7 picograms containing approximately two copies of the 3-billion-base-pair human genome—poses significant technical challenges for comprehensive genomic analysis [5–7]. To overcome this limitation, WGA has become essential for single-cell genomic studies [4, 7–10]. Various WGA techniques have been developed, each with distinct amplification mechanisms and characteristic error profiles. Multiple Displacement Amplification (MDA), introduced by Dean et al. [10], utilizes the highly processive phi29 DNA polymerase to achieve isothermal amplification with products exceeding 10 kb, though it suffers from pronounced amplification bias and chimera formation [11, 12]. Degenerate Oligonucleotide-Primed PCR (DOP-PCR), pioneered by Telenius et al. [13], employs thermocycling with degenerate primers to achieve more uniform coverage but generates shorter amplicons. Multiple Annealing and Looping-based Amplification Cycles (MALBAC) combines quasi-linear preamplification with exponential amplification to reduce bias [8], while Linear Amplification via Transposon Insertion (LIANTI) uses transposon insertion to create defined amplification origins, significantly improving uniformity and reducing artifacts [7]. More recently, Primary Template-directed Amplification (PTA) [14] and droplet-based MDA (dMDA) [15, 16] have emerged as promising alternatives that modify reaction conditions to suppress chimera formation, though these methods require specialized equipment and protocols that have limited their widespread adoption. These amplification methods can increase DNA content by several orders of magnitude (typically 1,000- to 10,000-fold), generating sufficient material for high-coverage sequencing necessary for reliable variant calling, copy number analysis, and SV detection [4, 17–21].

Accurate single-cell genomics is particularly critical for multiple applications where false-positive SVs can lead to incorrect biological conclusions. In cancer research, distinguishing genuine clonal evolution patterns from amplification artifacts is essential

for understanding tumor heterogeneity and therapeutic resistance [3]. In developmental biology, accurate detection of somatic mosaicism enables the reconstruction of lineage relationships and identification of pathogenic mutations in rare cell populations. For CRISPR-based genome editing, single-cell analysis with reliable SV detection is crucial for comprehensive assessment of off-target effects and ensuring genomic stability [14]. However, false-positive SVs introduced during amplification can confound these analyses, leading to misinterpretation of genomic rearrangements and their biological significance [4, 22].

Despite its critical role, WGA introduces systematic artifacts that significantly impact downstream analyses [7, 11, 12, 22, 23]. Chief among these are chimeric sequences—artificial DNA constructs formed through template switching during amplification. During MDA, the highly processive phi29 polymerase can dissociate from one genomic template and reinitiate synthesis on another, creating chimeric molecules that join DNA fragments from distant genomic loci into single amplified products [11]. These artifacts are particularly problematic for long-read sequencing technologies, where chimeric reads can span tens of kilobases and produce convincing—but false—evidence for genomic rearrangements.

Current computational tools to detect SVs from long-read data, including Sniffles2 [24, 25], DeBreak [26], SVIM [27], and cuteSV [28]. These methods typically employ read alignment analysis, split-read detection, and local assembly strategies to identify SV signatures [29]. However, distinguishing genuine biological SVs from WGA-induced chimeric artifacts remains challenging [23, 30–32].

Current computational approaches for identifying WGA-induced artifacts rely primarily on coverage-based metrics and read-pair orientation patterns [23, 30]. However, these heuristic methods often fail to distinguish genuine SVs from amplification artifacts, particularly when chimeric sequences exhibit complex rearrangement patterns, occur in repetitive genomic regions, or involve multiple genomic loci [31, 32]. This lack of robust, automated artifact detection has limited the reliability of SV analysis in single-cell studies and hindered the full realization of single-cell genomics’ potential for studying somatic mosaicism, tumor evolution, and rare cell populations.

The emergence of deep learning, particularly language models based on transformer architectures, has demonstrated remarkable success in genomics applications [33–36]. Recent genomic language models have shown the ability to learn complex sequence patterns and contextual relationships in DNA sequences, enabling improved performance in tasks such as regulatory element prediction, variant effect prediction, and functional annotation [36–38]. These models treat DNA sequences analogously to natural language, learning representations that capture both local motifs and long-range dependencies [33]. By training on large-scale genomic datasets, such models can internalize patterns of genuine biological sequences, including characteristic features of repetitive elements, chromatin structure, and sequence composition biases.

Here, we developed ChimeraLM, a genomic language model specifically designed to detect chimeric artifacts introduced by WGA. By leveraging deep learning to capture sequence patterns, structural features, and contextual information in genomic reads [33–36, 38], ChimeraLM effectively distinguishes genuine biological sequences from WGA-induced chimeric artifacts. We demonstrate that ChimeraLM achieves

139 superior performance compared to existing methods and substantially improves the
140 reliability of [SV](#) detection in single-cell genomic studies, thereby enabling accurate SV
141 analysis at single-cell resolution.

142 143 Results 144

145 Overview of ChimeraLM workflow and model architecture

146 Single-cell genomics relies on [WGA](#) to obtain sufficient DNA for sequencing (Fig. 1a).
147 The standard workflow includes single-cell isolation, DNA extraction, [WGA](#), long-read
148 sequencing (e.g., [Oxford Nanopore Technologies \(ONT\)](#)), base calling, and alignment
149 to the reference genome. During amplification, template-switching events introduce
150 artificial chimeric reads, resulting in alignment files that contain a mixture of authentic
151 and artificial sequences. In downstream analysis, these artifacts can mimic [SV](#) and
152 confound variant detection. To address this challenge, we developed ChimeraLM, a
153 [Genomic Language Model \(GLM\)](#) designed to integrate directly into this analysis
154 pipeline and distinguish biological reads from amplification-induced artifacts.

155 ChimeraLM functions as a pre-processing filter, operating after read alignment
156 but before [SV](#) detection. It evaluates each chimeric read—sequences with multiple
157 alignments to distant genomic locations—and classifies it as either biological (genuine)
158 or artificial ([WGA](#)-induced). This binary decision enables the retention of authentic
159 genomic sequences while removing amplification artifacts prior to variant calling. The
160 resulting high-confidence biological reads are then passed to conventional [SV](#) detection
161 algorithms for accurate identification of genomic rearrangements.

162 A high-confidence labeled dataset was required for supervised training of the model
163 (Fig. 1b; Extended Data Fig. 1a). We constructed this dataset using sequencing data
164 from the PC3 prostate cancer cell line, which provides both [WGA](#)-amplified and
165 non-amplified (bulk) genomic data. The key assumption is that bulk sequencing con-
166 tains only genuine genomic sequences, whereas [WGA](#) data includes both genuine and
167 artificial chimeras. Chimeric reads from the PC3 [WGA](#) PromethION dataset were sys-
168 tematically compared against three independent bulk datasets ([ONT](#) PromethION,
169 [ONT](#) MinION, and PacBio; see [Methods](#)). [WGA](#) reads whose chimeric structures were
170 absent from all three bulk datasets were labeled artificial. Conversely, [WGA](#) reads
171 with structures validated in one or more bulk datasets were labeled biological.

172 Application of this labeling strategy to the PC3 [WGA](#) data (Extended Data
173 Table 1) quantified the read distribution across these categories (Extended Data
174 Fig. 1b). We identified 12,670,396 chimeric reads with zero matches in the bulk ref-
175 erence, which were classified as artificial. Conversely, we identified a total of 293,180
176 reads validated as biological. This biological set was composed of reads matching one
177 (Match 1: 101,094 reads), two (Match 2: 190,309 reads), or all three (Match 3: 1,777
178 reads) of the bulk reference datasets. To construct a balanced training dataset, we
179 retained all 293,180 biological reads (combining Match 1, 2, and 3) and subsampled
180 an equal number (293,180) of artificial reads from the no-match category. This set was
181 augmented with 178,748 chimeric reads subsampled from the bulk datasets as posi-
182 tive controls. The final dataset of 765,108 labeled reads was partitioned into training
183 (70%), validation (20%), and internal test (10%) sets using stratified splitting.
184




Fig. 1 ChimeraLM workflow and architecture for detecting WGA artifacts in single-cell sequencing. (a) Single-cell genomic workflow and ChimeraLM integration. Single cells are isolated, followed by DNA extraction and WGA for genome amplification. WGA generates chimeric artifacts (red) through template switching during amplification, alongside biological reads (green). After nanopore sequencing, ChimeraLM classifies chimeric reads as biological or artificial, enabling downstream SV detection on clean reads. (b) Ground truth label generation for supervised learning. Chimeric reads from bulk sequencing data of the same cell line are compared against all chimeric reads from single-cell sequencing data. Reads that match bulk data are labeled as biological (green pathway), while non-matching reads are labeled as chimera artifacts (red pathway). This provides reliable training labels. (c) ChimeraLM architecture. Input DNA sequences (batch size B , sequence length L) are tokenized and encoded into hidden states $H \in \mathbb{R}^{L \times 256}$ through a backbone encoder (HyenaDNA [35]). Hyena operators capture long-range dependencies in genomic sequences. Attention pooling aggregates position-specific features using learned weights. Residual and multilayer perceptron (MLP) layers process pooled features, and a softmax layer outputs binary classification probabilities for biological versus artificial reads. (d) Attention pooling mechanism detail. The backbone encoder (HyenaDNA) transforms input sequences into hidden state $H \in \mathbb{R}^{L \times 256}$. Attention weights $\alpha \in \mathbb{R}^{L \times 1}$ are computed through linear layers, GELU activation, and softmax normalization, assigning importance scores to each nucleotide position. The weighted sum $h_{\text{pooled}} = \sum_{i=1}^L \alpha_i h_i$ produces the pooled output $h_{\text{pooled}} \in \mathbb{R}^{256}$, compressing variable-length sequences into fixed-dimensional representations. Created with BioRender.com.

The architecture of ChimeraLM (Fig. 1c) was specifically designed to learn from this dataset by operating directly on raw DNA sequences, bypassing conventional, feature-based classifiers. This design must address three primary technical challenges: (1) efficiently processing variable-length sequences of many kilobases, (2) simultaneously maintaining single-nucleotide resolution to detect the precise, abrupt compositional changes that define chimeric junctions, and (3) aggregating variable-length sequence representations into a consistent classification output.

ChimeraLM first addresses the need for high resolution by tokenizing input sequences at the single-nucleotide level. This base-pair precision is required to preserve the complete sequence information necessary for detecting chimeric junctions—the breakpoints where disparate genomic regions are artificially fused and which often exhibit abrupt compositional changes.

The architecture’s core employs Hyena operators [39], selected specifically to overcome the challenge of processing long DNA sequences. Traditional attention mechanisms scale quadratically with sequence length, making them computationally prohibitive for long-read data. Hyena operators, by contrast, achieve subquadratic scaling, enabling ChimeraLM to analyze full-length reads without fragmentation and thus preserve the structural context around chimeric junctions. To leverage existing genomic knowledge, we initialized the model with weights from HyenaDNA [35], a genomic foundation model pre-trained on diverse DNA sequences.

Finally, to produce a classification, the model employs an attention pooling mechanism to aggregate information across the entire variable-length read (Fig. 1d). This module computes learned, position-specific weights to identify which nucleotides—such as those at the junction boundary—are most informative for the classification decision. This weighted aggregation produces a fixed-dimensional representation, which is then processed through MLP components with residual connections. A final softmax layer outputs the probability scores for the biological versus artificial classes (see Methods). This end-to-end architecture enables ChimeraLM to learn directly from raw sequence data, discovering complex patterns that may not be apparent through rule-based algorithms.

ChimeraLM achieves high accuracy and reduces artifacts to near-bulk levels across platforms

We first evaluated ChimeraLM’s classification accuracy on the held-out test set (derived from the PromethION training data), which comprised reads with known biological or artificial status (Fig. 2a). The model achieved an F1 score of 0.81, reflecting balanced sensitivity and specificity in artifact detection. A recall of 0.95 indicates that 95% of true chimeric reads were correctly identified—critical for minimizing downstream false-positive structural variant calls—while a precision of 0.70 shows that the majority of reads flagged as chimeric were true artifacts. These results establish the model’s reliability for identifying amplification-induced artifacts in long-read data.

We next assessed its practical effectiveness on the full PC3 WGA datasets, comparing performance on the PromethION and MinION platforms (Fig. 2b). Bulk sequencing established a low baseline chimeric read rate (2.3% for PromethION; 2.5% for MinION). WGA dramatically increased this artifact load to 46.0% (PromethION)




Fig. 2 ChimeraLM accurately identifies and removes WGA-induced chimeric artifacts. (a) Classification performance on held-out test data. ChimeraLM achieves high recall (0.95) in identifying chimera artifacts while maintaining acceptable precision (0.70), yielding an F1 score of 0.81 for binary classification of biological versus artificial sequences. (b) Chimeric read reduction across sequencing platforms. Stacked bars show the proportion of chimeric (dark teal) and non-chimeric (light teal) reads in bulk sequencing, WGA-amplified samples, and ChimeraLM-filtered WGA samples. Data from PC3 cell line sequenced on PromethION (left) and MinION (right) platforms demonstrate that ChimeraLM reduces chimeric read frequencies from 46.0% to 4.9% (PromethION) and from 23.0% to 1.5% (MinION), approaching bulk levels (2.3% and 2.5%, respectively). (c,d) Benchmarking against existing methods. ChimeraLM achieves approximately 90% reduction in chimeric reads on both PromethION (c) and MinION (d) platforms, whereas existing computational tools SACRA and 3rd-ChimeraMiner show no detectable reduction in chimeric content.

and 23.0% (MinION). After ChimeraLM filtering, chimeric content dropped to 4.9% on PromethION and 1.5% on MinION—representing 10- to 15-fold reductions—while retaining 15.8 million and 5.6 million biological reads. This restoration to near-bulk quality demonstrates that ChimeraLM effectively separates genuine genomic reads from WGA-induced artifacts.

We then benchmarked ChimeraLM against existing computational tools for detecting amplification-induced chimeras, SACRA [30] and 3rd-ChimeraMiner [23] (Fig. 2c,d). When applied to the same PromethION and MinION WGA data, ChimeraLM achieved an approximately 90% reduction in chimeric reads on both platforms. In stark contrast, neither SACRA nor 3rd-ChimeraMiner showed any detectable reduction in chimeric content (0% reduction).

Together, these results demonstrate a robust and generalizable performance. The strong filtering on the MinION dataset (Fig. 2b) is particularly noteworthy, as this dataset served as a completely independent test set; the model was trained exclusively on PromethION data. This cross-platform generalization, combined with the

high recall on the internal test set (Fig. 2a) and the clear superiority over existing tools (Fig. 2c,d), indicates that ChimeraLM learns fundamental, generalizable sequence features of WGA-induced artifacts rather than platform-specific signatures.




Fig. 3 ChimeraLM improves structural variant detection accuracy. (a) Construction of high-confidence SV reference dataset. PC3 bulk DNA was sequenced on multiple platforms (ONT PromethION and MinION, Illumina HiSeq) and analyzed with multiple SV calling algorithms. SV events detected by ≥ 2 callers on the same platform were retained. Events supported by both long-read and short-read platforms were designated as high-confidence gold standard SVs. (b,c) SV validation against multi-platform gold standard. Stacked bars show total SV calls (log scale, numbers above bars) classified as gold standard-supported (dark teal) or unsupported (light teal) for PromethION (b) and MinION (c). ChimeraLM substantially reduces unsupported SV calls while preserving gold standard events. (d,e) SV validation against long-read bulk sequencing (ONT PromethION and MinION). Stacked bars show SV calls classified as bulk-supported (dark teal) or unsupported (light teal) for PromethION (d) and MinION (e). Long-read bulk data from the same platform provides platform-matched validation, capturing true variants that may be specific to long-read detection. (f,g) SV type distribution across processing methods. Bar charts show the number of detected SVs by type: deletion (DEL) (green), duplication (DUP) (orange), insertion (INS) (blue), inversion (INV) (pink), and translocation (TRA) (light green). Unfiltered WGA data shows elevated counts across all types, particularly INVs and TRAs, which are reduced to bulk-like levels after ChimeraLM filtering. (h,i) Composition of chimeric artifact-supported SVs. Pie charts show the proportion of SV types among events supported exclusively by reads classified as chimeric artifacts in unfiltered WGA data for PromethION (h) and MinION (i). These represent false-positive SV calls that would be eliminated by ChimeraLM.

ChimeraLM substantially reduces false-positive structural variant calls

Accurate SV detection is essential for understanding genomic diversity and disease mechanisms in single cells. However, WGA-induced chimeric artifacts can be misidentified as genuine SVs, leading to incorrect biological conclusions. To quantify ChimeraLM's impact on SV calling accuracy, we compared variant calls from unfiltered WGA data and ChimeraLM-filtered data against two independent reference standards (Fig. 3).

We first established a high-confidence gold standard SV dataset by integrating results from bulk PC3 DNA sequenced on multiple platforms (ONT PromethION, ONT MinION, and Illumina HiSeq) and analyzed with multiple SV callers (Fig. 3a; Extended Data Table 1). SVs detected by ≥ 2 callers on the same platform and supported by both long-read and short-read data were retained as gold-standard events, ensuring high specificity across technologies.

Comparison against this gold standard revealed that unfiltered WGA data contained extensive false-positive SVs (Fig. 3b,c). On PromethION, raw WGA data produced 3.6 million SV calls, of which only 8,815 (0.24%) matched gold standard events—indicating that over 99% were artifacts. After ChimeraLM filtering, total calls dropped to 305,570 while retaining 8,067 true events, raising the validation rate to 2.64% (11-fold improvement) and preserving 91.5% of true variants. MinION data showed similar results, with calls reduced from 451,390 to 38,432 and the validation rate increasing from 1.59% to 16.3% (10-fold improvement) while retaining 87.2% of true variants. These results highlight ChimeraLM's ability to remove spurious SV calls while maintaining biological sensitivity.

To complement this stringent validation, we next performed platform-matched bulk validation, comparing WGA-derived SV calls against long-read bulk sequencing from the same platform (Fig. 3d,e). This reference captures true SVs that may be missed by short-read data, providing a more inclusive measure of recall. Under this benchmark, ChimeraLM increased validation rates from 0.38% to 3.24% on PromethION (8.5-fold improvement) and from 1.25% to 12.79% on MinION (10-fold improvement), while retaining 71.9% and 87.1% of bulk-supported events, respectively. The consistent improvements across independent datasets demonstrate that ChimeraLM effectively suppresses WGA-induced artifacts without sacrificing detection of genuine SVs.

Together, these analyses demonstrate that ChimeraLM reduces false-positive SV calls by 8–16 fold while preserving 72–92% of true variants, resulting in a substantial enhancement of the signal-to-noise ratio in single-cell SV discovery. By restoring near-bulk specificity and maintaining robust sensitivity, ChimeraLM enables more accurate and interpretable downstream genomic analyses.

ChimeraLM restores unbiased SV-type distributions and characterizes artifact composition

Amplification artifacts can distort the apparent spectrum of SVs, often inflating specific SV types. To evaluate whether ChimeraLM effectively corrects such distortions,

we compared **SV** type distributions across bulk, unfiltered **WGA**, and ChimeraLM-filtered datasets (Fig. 3f,g). Bulk sequencing showed relatively balanced proportions of **DELs**, **DUPs**, **INSSs**, **INVs**, and **TRAs**. In contrast, unfiltered **WGA** data exhibited a dramatic overrepresentation of **INVs** on both PromethION and MinION platforms, consistent with pervasive amplification artifacts. After ChimeraLM filtering, these distributions were largely restored toward bulk-like profiles: excessive **INVs** were markedly reduced while other **SV** categories remained stable. This shift reflects selective removal of artifact-supported **INVs** rather than indiscriminate loss of genuine inversion signals, demonstrating high specificity in distinguishing chimeric from biological reads.

To investigate the basis of this normalization, we analyzed **SV** calls supported exclusively by reads classified as chimeric by ChimeraLM (Fig. 3h,i). These artifact-supported events were overwhelmingly dominated by **INVs**, comprising 88.4% on PromethION and 92.4% on MinION. This pattern is consistent with template-switching junctions that produce inversion-like alignment signatures. Smaller fractions of **DELs** (5.1% and 3.8%), **DUPs** (3.4% and 2.4%), and **INSSs** (3.0% and 1.4%) were also observed, demonstrating that **WGA**-induced chimeras can mimic diverse **SV** categories rather than only **INVs**.

This characterization has important implications for single-cell genomics. Although **INVs** are the predominant artifact type, the coexistence of **DELs**, **DUPs**, and **INSSs** among chimeric events indicates that comprehensive filtering—rather than inversion-specific correction—is essential for accurate **SV** detection. Without ChimeraLM filtering, single-cell **SV** analyses would be confounded not only by false-positive **INVs** but also by other artifact-associated variants [31, 32]. By restoring biologically representative **SV** type distributions, ChimeraLM enables robust and interpretable characterization of structural variation in single cells without distortion from **WGA**-induced artifacts.




Fig. 4 ChimeraLM attention weights can localize to chimeric junction regions.

(a,d) Attention weight profiles for two representative chimeric reads. Upper panels show attention weights per sequence position (blue line) and mean attention (dashed line). Red vertical lines mark chimeric junction positions, with pink shading indicating junction region (± 50 bp). Purple arrows show reported SV breakpoints. Lower panels illustrate read alignments: reads (top bars) show orientation transitions at junctions (green = forward, red = reverse-complemented, arrows indicate strand), while reference genome (bottom bars) maintains continuous forward orientation. Gray regions connect aligned segments. (b,e) Quantitative attention analysis. Box plots show significantly elevated attention weights in junction region versus non-junction regions for both examples ($p = 5.3 \times 10^{-14}$ and $p = 6.8 \times 10^{-15}$, respectively; Wilcoxon rank-sum test). (c,f) Proposed chimera formation mechanisms. Step 1: Original DNA fragments from distant genomic loci exist in forward orientation. Step 2: During WGA, one or both fragments may undergo random reverse-complementation. Step 3: Template switching joins the fragments with discordant orientations, creating chimeric artifacts. The two examples illustrate different orientation patterns (forward-to-reverse vs reverse-to-forward transitions) arising from random strand selection during amplification.

ChimeraLM provides interpretable classification through attention visualization

We next investigated whether ChimeraLM’s attention mechanism highlights biologically meaningful regions within sequencing reads (Fig. 4).

For representative chimeric reads, attention weight profiles showed low baseline values across most positions but pronounced peaks at junction regions where template switching artificially joins DNA fragments from distinct genomic loci (Fig. 4a,d). These peaks coincided precisely with alignment breakpoints characterized by orientation changes between adjacent read segments—the defining signature of WGA-induced chimeric artifacts.

Quantitative analysis confirmed that attention weights within junction regions (± 50 bp) were significantly higher than those in non-junction regions (Wilcoxon rank-sum test, $p = 5.3 \times 10^{-14}$ and $p = 6.8 \times 10^{-15}$) (Fig. 4b,e). Such localization indicates that ChimeraLM learns mechanistically relevant features associated with artificial junction formation rather than relying on spurious correlations.

Schematic reconstruction of the amplification process further supports this interpretation (Fig. 4c,f). During WGA, DNA fragments from distant genomic loci may undergo random strand orientation changes before being joined by template switching. This process produces artificial junctions with discordant orientations—forward-to-reverse or reverse-to-forward—that generate inversion-like alignment signatures and are effectively recognized by the model’s attention peaks.

Together, these analyses demonstrate that ChimeraLM’s attention mechanism can localize chimeric junctions at single-base resolution and capture the underlying orientation discontinuities that define WGA-induced artifacts.

Discussion

WGA has enabled genomic analysis from single cells but introduces chimeric artifacts that compromise SV detection. ChimeraLM addresses this challenge through sequence-level classification of biological versus artificial reads, substantially improving SV calling accuracy before downstream analysis. This upstream filtering strategy—removing problematic sequences at the read level rather than correcting errors post hoc—provides a practical solution for single-cell genomics laboratories.

Our results demonstrate several key advantages of ChimeraLM for long-read single-cell sequencing. The method achieves approximately 90% reduction in chimeric reads across nanopore platforms while retaining 87–92% of true SVs. It reduces false-positive SV calls by 8–16 fold, enabling researchers to focus on biologically relevant variants without manually filtering thousands of artifacts. Moreover, ChimeraLM performs consistently across PromethION and MinION without platform-specific retraining, indicating that it captures generalizable sequence features of WGA-induced chimeras. These results underscore the model’s robustness across diverse datasets and sequencing conditions.

ChimeraLM’s effectiveness reflects the ability of deep learning models to capture complex sequence patterns that are difficult to encode in rule-based filters. Traditional quality control methods rely on predefined metrics such as mapping quality or read

depth [23, 30], which may not effectively distinguish chimeric artifacts from biological reads. By learning directly from data, ChimeraLM discovers subtle compositional and structural features that differentiate authentic genomic sequences from amplification artifacts. Furthermore, the model offers interpretability through attention visualization, allowing researchers to examine which sequence regions drive classification. Attention weights can concentrate sharply at junctions where template switching joins DNA fragments from distinct loci, matching the known mechanism of chimera formation. Some reads show more diffuse attention distributions, suggesting that ChimeraLM integrates multiple complementary cues—such as junction orientation, compositional biases, and local sequence context—to classify diverse artifact types. This interpretability builds confidence in the model’s predictions and provides a lens for probing the molecular processes underlying amplification-induced artifacts.

The improved reliability of SV detection has direct implications for single-cell genomics. Studies of chromosomal instability, clonal evolution, and SV burden in individual cells have long been constrained by high false-positive rates in WGA data [31, 32]. ChimeraLM enables more confident identification of genuine SVs, supporting research in cancer genomics, developmental biology, and aging where single-cell resolution is essential for understanding cellular heterogeneity. Although the current model processes reads independently, integrating additional contextual features—such as coverage, mate-pair, or phasing information—could further enhance accuracy. Graphics Processing Unit (GPU) resources are recommended for large-scale datasets, while Central Processing Unit (CPU) inference remains feasible for smaller studies; runtime optimization and model compression may improve accessibility for broader use.

Future work should prioritize validation across diverse biological contexts in long-read single-cell sequencing. Testing on multiple cell types (primary, stem, or immune cells) and WGA protocols will establish generalizability. The interpretability of attention-based models could also be leveraged to investigate mechanisms of chimera formation: large-scale analysis of attention patterns may reveal recurrent sequence motifs or genomic contexts associated with template switching, guiding the development of improved amplification protocols. More broadly, ChimeraLM illustrates the potential of GLMs for data quality control applications [35]. Architectural innovations such as the Hyena operator for efficient long-range modeling [39] may have utility beyond chimera detection, addressing challenges such as contamination, adapter artifacts, and systematic sequencing errors.

Looking ahead, ChimeraLM’s framework could extend beyond single-cell genomics to address quality control challenges in other amplification-dependent technologies, including cell-free DNA analysis, ancient DNA studies, and metagenomic sequencing from low-biomass samples. The model’s interpretability through attention visualization also opens opportunities for mechanistic studies of polymerase fidelity and template-switching dynamics across different amplification protocols. Furthermore, integration with emerging single-cell multi-omics platforms could enable simultaneous quality control across genomic, transcriptomic, and epigenomic data layers, providing a unified framework for artifact detection in complex single-cell experiments.

ChimeraLM thus provides a practical and interpretable framework for improving long-read single-cell genomic data quality. By removing WGA-induced chimeric artifacts at the read level and revealing the mechanistic features that drive them, the method not only enhances SV detection reliability but also deepens understanding of amplification-induced bias in single-cell genomics.

Methods

Cell culture, single-clone preparation, and nanopore sequencing

Cell culture and single-clone establishment

PC3 prostate cancer cells (ATCC[®] CRL-1435[™]) were cultured in RPMI-1640 medium supplemented with 10% fetal bovine serum and 1% penicillin–streptomycin at 37 °C with 5% CO₂. To minimize biological heterogeneity, a monoclonal population was established by serial dilution in 96-well plates, ensuring that each culture originated from a single cell. Mycoplasma contamination was routinely tested and confirmed negative prior to DNA extraction.

DNA extraction and whole-genome amplification

From the monoclonal population, two types of DNA samples were prepared: a bulk (non-amplified) control and ten single-cell MDA-amplified genomes. Bulk high-molecular-weight DNA was extracted using the Monarch[®] HMW DNA Extraction Kit for Cells & Blood (New England Biolabs). Individual cells were isolated using 1CellDish-60 mm (iBiochips) and amplified using the REPLI-g Advanced DNA Single Cell Kit (Qiagen) following the manufacturer’s protocol. DNA concentration and fragment integrity were assessed with a Qubit 4 fluorometer and Agilent TapeStation (DNA 1000/5000 ScreenTape). Only samples meeting quality standards were used for library construction.

Nanopore library preparation and sequencing

Sequencing libraries were prepared using the ONT Ligation Sequencing Kit V14 (SQK-LSK114) and sequenced on MinION Mk1C or PromethION P2 Solo devices with R10.4.1 flow cells according to the manufacturer’s genomic DNA workflow. Because all single-cell samples originated from the same monoclonal lineage, observed differences between amplified and bulk data primarily reflect MDA-induced artifacts rather than biological variation, providing a controlled experimental setting for downstream analyses.

Basecalling and read processing

Raw signal files (POD5) were basecalled using Dorado v0.5.0 with the high-accuracy model dna_r10.4.1_e8.2.400bps_hac@v4.3.0 [40]. Reads with mean quality < 10 or length < 500 bp were removed. Residual adapters and concatemers were trimmed using Cutadapt v4.0 [41] in two-pass error-tolerant mode. Cleaned reads were aligned to the GRCh38.p13 reference genome using minimap2 v2.26 (map-ont preset) [42]. Resulting BAM files were sorted and indexed with SAMtools v1.16 [43]. Read length

691 and mapping statistics were calculated using NanoPlot v1.46.1 [44]. All samples were
692 processed under identical parameters to ensure consistency across datasets.

693
694 *Chimeric read identification*

695 Chimeric reads were identified based on the presence of supplementary alignments in
696 BAM files using the [Supplementary Alignment \(SA\)](#) tag. The [SA](#) tag indicates that
697 a read has additional alignments beyond the primary alignment, which is character-
698 istic of chimeric sequences that map to multiple distant genomic locations. To ensure
699 accurate identification, we applied stringent filtering criteria: reads were classified as
700 chimeric only if they (1) were not unmapped, (2) contained the [SA](#) tag, (3) were not
701 secondary alignments, and (4) were not supplementary alignments themselves. This
702 filtering approach ensures that only primary alignments with supplementary mapping
703 evidence are considered chimeric, avoiding double-counting of the same chimeric event
704 and excluding low-quality or ambiguous alignments. Reads without the [SA](#) tag (single
705 continuous alignments) were classified as non-chimeric. This approach leverages the
706 standard BAM format specification to reliably identify reads with complex alignment
707 patterns.

708
709 **Training data construction**

710
711 *Data generation and sources*

712 To construct the training dataset, we generated [WGA](#) and bulk sequencing data from
713 PC3 cells. The [WGA](#) sample was amplified and sequenced on the PromethION P2 plat-
714 form ([ONT](#)), while three independent bulk datasets were produced from non-amplified
715 genomic DNA: bulk PromethION P2, bulk MinION Mk1c ([ONT](#)), and bulk PacBio.
716 These bulk datasets represent authentic biological sequences free from amplification-
717 induced artifacts. In contrast, [WGA](#) sequencing includes both genuine genomic reads
718 and artificial chimeras introduced during the amplification process. An additional
719 [WGA](#) dataset sequenced on the MinION Mk1c platform was reserved exclusively as
720 an independent test set for cross-platform evaluation.

721
722 *Ground truth annotation and class definition*

723 Ground truth labels were established by systematically comparing chimeric reads from
724 the [WGA](#) PromethION P2 dataset against those from the three bulk datasets. For
725 each [WGA](#) chimeric read, all alignment segments—defined by their genomic start
726 and end coordinates—were compared to the corresponding segments of bulk chimeric
727 reads. A [WGA](#) read was labeled as biological if every segment matched at least one
728 bulk chimeric read within a 1 kb positional tolerance, indicating that the structural
729 configuration is also present in non-amplified DNA. Reads lacking any matching pat-
730 tern across all bulk datasets were labeled as artificial chimeras, presumed to arise
731 from the amplification process. To ensure balanced class representation, additional
732 chimeric reads were randomly sampled from the bulk datasets and labeled as biologi-
733 cal, as these reads originate from genuine genomic rearrangements such as true [SVs](#).
734 The final labeled dataset combined the annotated [WGA](#) PromethION P2 reads with
735

736

the subsampled bulk chimeric reads and was subsequently partitioned into training, validation, and test sets as described below.

Dataset partitioning and cross-platform validation

The combined labeled dataset, derived from WGA PromethION P2 and bulk sequencing data, was divided into training (70%), validation (20%), and internal test (10%) sets using stratified random sampling to maintain class balance. These subsets were used respectively for model training, hyperparameter tuning, and performance evaluation on data from the same sequencing platform.

To evaluate cross-platform generalization, the complete WGA MinION Mk1c dataset was reserved as an independent external test set. This dataset, generated on a different nanopore platform, was never used during model training or internal testing. This two-level evaluation design allowed us to test whether ChimeraLM captures general sequence features of amplification-induced chimeras rather than platform-specific artifacts.

Model architecture

Backbone encoder

ChimeraLM employs the pre-trained HyenaDNA model [35] as its backbone encoder. This model was pre-trained on large-scale genomic data and provides robust sequence representations. DNA sequences are tokenized at single-nucleotide resolution, with each base (A, C, G, T, N) mapped to a unique integer token (7, 8, 9, 10, 11, respectively). Special tokens include [CLS]=0, [PAD]=4, and others for sequence processing. Input sequences are truncated at 32,768 bp or padded to enable batch processing.

For a tokenized input sequence $\mathbf{x} \in \mathbb{Z}^L$, the HyenaDNA backbone generates contextualized hidden representations:

$$\mathbf{H} = \text{HyenaDNA}(\mathbf{x}) \in \mathbb{R}^{L \times 256}$$

where $\mathbf{H} = (\mathbf{h}_1, \mathbf{h}_2, \dots, \mathbf{h}_L)$ represents position-wise hidden states with dimension 256. The Hyena operators [39] efficiently capture both local sequence motifs and long-range dependencies essential for distinguishing biological sequences from chimeric artifacts.

Attention pooling

To aggregate variable-length sequence representations into fixed-size vectors, ChimeraLM implements attention-based pooling. For hidden states $\mathbf{H} \in \mathbb{R}^{L \times 256}$, attention weights are computed through a two-layer network:

$$\mathbf{e} = \text{GELU}(\text{Linear}_{256 \rightarrow 256}(\mathbf{H})) \in \mathbb{R}^{L \times 256}$$

$$\mathbf{s} = \text{Linear}_{256 \rightarrow 1}(\mathbf{e}) \in \mathbb{R}^{L \times 1}$$

$$\boldsymbol{\alpha} = \text{softmax}(\mathbf{s}) \in \mathbb{R}^{L \times 1}$$

783 The pooled representation is the weighted sum of hidden states:

784

785

786

787

$$\mathbf{h}_{\text{pooled}} = \sum_{i=1}^L \alpha_i \mathbf{h}_i \in \mathbb{R}^{256}$$

788

789

790

791

This mechanism assigns learned importance weights to each sequence position, enabling the model to focus on informative regions while accommodating natural variability in read lengths.

792

793

Classification head

794

795

The pooled representation is processed through a [MLP](#) with residual connections. The first layer expands dimensionality:

796

797

$$\mathbf{f}_1 = \text{Dropout}_{0.1}(\text{GELU}(\text{Linear}_{256 \rightarrow 512}(\mathbf{h}_{\text{pooled}}))) \in \mathbb{R}^{512}$$

798

799

Subsequent residual blocks with input $\mathbf{f}_{\text{in}} \in \mathbb{R}^{512}$ compute:

800

801

$$\mathbf{f}_{\text{out}} = \text{Dropout}_{0.1}(\text{Linear}_{512 \rightarrow 512}(\text{GELU}(\text{Linear}_{512 \rightarrow 512}(\mathbf{f}_{\text{in}})))) + \mathbf{f}_{\text{in}}$$

802

803

804

where the skip connection enables stable gradient flow during training. The final layer produces binary classification logits:

805

806

$$\mathbf{z} = [z_0, z_1] = \text{Linear}_{512 \rightarrow 2}(\mathbf{f}_{\text{final}}) \in \mathbb{R}^2$$

807

808

809

where z_0 and z_1 represent logits for biological and artificial chimeric classes, respectively. During inference, the predicted class is $\hat{y} = \text{argmax}_{i \in \{0,1\}} z_i$.

810

811

Model summary

812

813

814

815

816

817

The complete ChimeraLM pipeline processes DNA sequences through: (1) single-nucleotide tokenization, (2) HyenaDNA backbone encoding to generate contextualized representations, (3) attention pooling to aggregate position-specific features, (4) [MLP](#) layers with residual connections to learn classification features, and (5) binary classification output. The entire model is trained end-to-end using labeled [WGA](#) and bulk sequencing data.

818

819

820

821

Training configuration

822

823

824

825

826

827

828

ChimeraLM was trained using PyTorch [\[45\]](#) and PyTorch Lightning [\[46\]](#) frameworks. Input sequences were tokenized using the tokenizer with maximum sequence length of 32,768 bp. Sequences longer than this threshold were truncated; shorter sequences were padded to enable batch processing. Training employed mixed-precision computation (bf16) to accelerate training while maintaining numerical stability.

Optimization procedure

We used the AdamW optimizer [47] with learning rate $\eta = 1 \times 10^{-4}$ and weight decay $\lambda = 0.01$. AdamW implements adaptive learning rates with decoupled weight decay, combining the benefits of Adam optimization with proper L2 regularization. A ReduceLROnPlateau scheduler dynamically adjusted the learning rate based on validation loss, reducing it by a factor of 0.1 when no improvement occurred for 10 consecutive epochs. Early stopping with patience of 10 epochs prevented overfitting by terminating training when validation performance plateaued. A fixed random seed (12345) ensured reproducibility across training runs.

The training objective used cross-entropy loss for binary classification. For a training example with true class label $y \in \{0, 1\}$ and model logits $\mathbf{z} = [z_0, z_1]$, the loss is:

$$\mathcal{L}(\mathbf{z}, y) = -\log \left(\frac{\exp(z_y)}{\exp(z_0) + \exp(z_1)} \right) = -z_y + \log(\exp(z_0) + \exp(z_1))$$

where z_0 and z_1 represent logits for biological and artificial chimeric classes, respectively.

Training implementation

Training used batch size of 16 sequences with 30 parallel data loading workers. GPU acceleration was employed for efficient processing, with training typically requiring 96-120 hours depending on dataset size. Model checkpointing saved the best-performing model based on validation metrics. Configuration management used Hydra [48] to enable reproducible experimentation.

Model evaluation

Performance was monitored using accuracy, precision, recall, and F1 score on the validation set after each epoch:

$$\begin{aligned} \text{Precision} &= \frac{\text{TP}}{\text{TP} + \text{FP}}, & \text{Recall} &= \frac{\text{TP}}{\text{TP} + \text{FN}} \\ \text{F1} &= \frac{2 \times \text{Precision} \times \text{Recall}}{\text{Precision} + \text{Recall}}, & \text{Accuracy} &= \frac{\text{TP} + \text{TN}}{\text{TP} + \text{TN} + \text{FP} + \text{FN}} \end{aligned}$$

where TP (true positives) are chimeric reads correctly classified as artificial, TN (true negatives) are biological reads correctly classified as biological, FP (false positives) are biological reads misclassified as artificial, and FN (false negatives) are chimeric reads misclassified as biological. Final model selection was based on best validation performance as determined by early stopping.

Model inference and application

Inference pipeline

To apply ChimeraLM to new WGA sequencing data, the model takes a BAM file as input. Chimeric reads are identified using SA tags and filtered to exclude unmapped, secondary, or supplementary alignments. Each chimeric read sequence is tokenized

875 using the tokenizer (maximum length 32,768 bp, with truncation or padding as
876 needed). The trained model processes sequences in batches, generating two logits
877 $[z_0, z_1]$ for each read corresponding to biological and artificial chimeric classes. Clas-
878 sification is determined by $\hat{y} = \text{argmax}(z_0, z_1)$. ChimeraLM outputs a filtered BAM
879 file containing only reads classified as biological, which can be directly used for
880 downstream analyses including SV calling.

881

882 Performance evaluation

883

884 *Test set evaluation*

885 Final model performance was evaluated on the held-out test set and the independent
886 MinION Mk1c dataset. Metrics (precision, recall, F1 score, accuracy) were computed
887 as described in the training section, where true positives represent chimeric reads
888 correctly classified as artificial and true negatives represent biological reads correctly
889 classified as biological.

890

891 *SV calling*

892 SVs were called using multiple tools to ensure comprehensive detection. For long-
893 read data (ONT PromethION P2 and MinION Mk1c), we used Sniffles v2.5 [24, 25],
894 DeBreak v1.2 [26], SVIM v2.0.0 [27], and cuteSV v2.1.1 [28]. For short-read data of the
895 PC3 cell line, we used both the CCLE Illumina whole-genome sequencing dataset and
896 the PRJNA361315 Illumina WGS dataset, processed with Manta v1.6.0 [49], DELLY
897 v1.5.0 [50], and SvABA v1.1.0 [51]. All tools were executed with default recommended
898 parameters.

899

900 *Gold standard SV dataset construction*

901 A high-confidence gold standard SV dataset was generated from bulk PC3 sequencing
902 data to evaluate the impact of ChimeraLM on SV detection accuracy (Fig. 3a). All
903 SV comparison and breakpoint correction were performed using OctopusSV v0.2.3 [52].
904 We used four datasets: bulk MinION Mk1c, bulk PromethION P2, the CCLE Illumina
905 WGS dataset, and the PRJNA361315 Illumina WGS dataset. Within each dataset, SV
906 events supported by at least two independent callers were retained. Variants supported
907 by two or more datasets were designated as gold standard SVs for benchmarking.

908

909 *SV benchmarking analysis*

910 To assess the impact of ChimeraLM on SV calling accuracy, we compared SV calls from
911 unfiltered WGA data and ChimeraLM-filtered WGA data against two references: (1)
912 the stringent multi-platform gold standard dataset, and (2) platform-matched long-
913 read bulk sequencing data. Benchmarking was performed using Truvari v4.2.2 [53]
914 with default parameters. SVs were considered supported if they matched reference
915 variants within the defined breakpoint tolerance. Validation rates were calculated as
916 the proportion of called SVs supported by the reference. This dual benchmarking
917 strategy quantifies both improvements in detecting high-confidence multi-platform
918 SVs and the retention of platform-specific true variants.


919

920

Benchmarking against existing methods	921
ChimeraLM was compared to two existing computational methods for detecting	922
amplification-induced chimeric artifacts: SACRA [30] (GitHub commit 9a2607e) and	923
3rd-ChimeraMiner [23] (GitHub commit 04b5233). Both tools were applied to WGA	924
data from PromethION P2 and MinION Mk1c platforms using default parameters as	925
recommended in their documentation. Performance was evaluated by measuring the	926
percentage reduction in chimeric reads relative to unprocessed WGA data. Chimeric	927
reads were identified using WGA tag-based alignment criteria (reads with SA tags	928
indicating split alignments), and reduction rates were calculated as the proportion of	929
chimeric reads removed by each method.	930
	931
Attention weight analysis	932
To investigate ChimeraLM’s interpretability, we analyzed attention weights from	933
the pooling mechanism for representative chimeric reads. Attention weights indicate	934
the relative importance assigned to each sequence position during classification. For	935
selected reads, we extracted per-position attention weights and visualized them along-	936
side read alignments to identify whether the model focuses on mechanistically relevant	937
regions.	938
Chimeric junction positions were identified from alignment data (defined by break-	939
points in SA tags). A window of ± 50 bp surrounding each junction was designated as	940
the junction region. Attention weights within junction region were compared to non-	941
junction regions using the Wilcoxon rank-sum test [54], with statistical significance	942
assessed at $p < 0.001$.	943
	944
Data visualization	945
Figures were generated using Python with Matplotlib [55] and Seaborn [56].	946
	947
Computing resources	948
Computations were performed on a High Performance Computing (HPC) server with	949
64-core Intel Xeon Gold 6338 CPU, 256 GB RAM, and two NVIDIA A100 GPUs (80	950
GB memory each).	951
	952
Supplementary information.	953
	954
Acknowledgements. We thank Tingyou Wang for guidance on figure preparation.	955
This project was supported in part by NIH grants R35GM142441 and R01CA259388	956
awarded to RY.	957
	958
Declarations	959
	960
Author Contributions. YL, QG and RY designed the study. YL and QG per-	961
formed the analysis. QG performed the experiments. YL and QG designed and	962
implemented the model. YL built the command-line tool and documentation. YL, QG	963
and RY wrote the manuscript. RY supervised this work.	964
	965
	966

Extended Data Table 1 Sequencing and alignment statistics of PC3

Sample	Platform	Reads ($\times 10^6$)	Total bases (Gb)	Total bases aligned (Gb)	Fraction aligned	Mean length (bp)	Mean quality (Q)	Average identity (%)
WGA	MinION	9.11	14.6	10.4	0.7	1,603	14.3	97.6
WGA	PromethION	44.69	128.2	69.2	0.5	2,869	14.5	96.1
Bulk	MinION	0.97	8.1	7.1	0.9	8,310	17.2	97.3
Bulk	PromethION	8.00	69.9	62.4	0.9	8,732	18.5	97.7



Extended Data Fig. 1 Training dataset construction and ground-truth labeling strategy for PC3 cell line. (a) Schematic workflow for generating labeled training data. WGA PromethION data containing both biological and artificial chimeric reads is compared against three independent bulk sequencing datasets from the same cell line (PromethION, MinION, and PacBio platforms). Chimeric reads are classified through systematic matching: reads with no matches across all bulk datasets (Match 0) are labeled as artificial chimeras (positive class, red); reads matching one or more bulk datasets (Match 1–3) are labeled as biological reads (negative class, green), along with chimeric reads sampled directly from bulk data. The combined labeled dataset undergoes stratified random sampling to generate training (70%), validation (20%), and testing (10%) sets for model development. The WGA MinION dataset is reserved as an independent cross-platform generalization test set. (b) Distribution of chimeric read matches between WGA and bulk sequencing datasets. Bar chart showing the number of chimeric reads (y-axis, log scale) grouped by how many bulk datasets (x-axis) contained matching chimeric structures when comparing WGA PromethION reads against bulk sequencing data. “Match 0” indicates reads with no matches in any bulk dataset (classified as artificial chimeras, $\sim 10^7$ reads), whereas “Match 1–3” indicate reads with matches in one, two, or all three bulk datasets (classified as biological reads, $\sim 10^5$ reads each). Color-coded boxes below bars indicate which bulk platforms validated each read category: PromethION (light blue), MinION (white), and PacBio (white); green boxes indicate platform-specific validation. The substantial imbalance between Match 0 ($\sim 10^7$) and Match 1–3 categories ($\sim 10^5$ each) reflects the high prevalence of WGA-induced artifacts, necessitating balanced subsampling for supervised learning.

Data Availability. The raw sequencing data generated in this study have been deposited in the NCBI Sequence Read Archive (SRA) under BioProject accession

PRJNA1354861. The dataset includes Oxford Nanopore long-read whole-genome sequencing of PC3 prostate cancer cells and MDA-amplified single-cell derivatives. The individual SRA accessions are as follows: PC3 bulk (MinION Mk1C), SRR35904028; PC3 bulk (PromethION P2), SRR35904029; PC3 10-cell WGA (MinION Mk1C), SRR35904026; PC3 10-cell WGA (PromethION P2), SRR35904027. We can access the data at the following link: <https://dataview.ncbi.nlm.nih.gov/object/PRJNA1354861?reviewer=viej6cv6mgbli3n7a9a5k1bsb3>

Code Availability. ChimeraLM, implemented in Python, is open source and available on GitHub (<https://github.com/ylab-hi/ChimeraLM>) under the Apache License, Version 2.0. The package can be installed via PyPI (<https://pypi.org/project/chimeralm>) using pip, with wheel distributions provided for Windows, Linux, and macOS to ensure easy cross-platform installation. An interactive demo is available on Hugging Face (<https://huggingface.co/spaces/yangliz5/ChimeraLM>), allowing users to test DeepChopper’s functionality without local installation. For large-scale analyses, we recommend using ChimeraLM on systems with GPU acceleration. Detailed system requirements and optimization guidelines are available in the repository’s documentation (<https://ylab-hi.github.io/ChimeraLM/>).

Conflict of interest. RY has served as an advisor/consultant for Tempus AI, Inc. This relationship is unrelated to and did not influence the research presented in this study.

Acronyms

CPU Central Processing Unit 14

DEL deletion 9, 11

dMDA droplet-based MDA 2

DOP-PCR Degenerate Oligonucleotide-Primed PCR 2

DUP duplication 9, 11

GLM Genomic Language Model 4, 14

GPU Graphics Processing Unit 14, 19, 21, 23

HPC High Performance Computing 21

INS insertion 9, 11

INV inversion 1, 9, 11

LIANTI Linear Amplification via Transposon Insertion 2

MALBAC Multiple Annealing and Looping-based Amplification Cycles 2

MDA Multiple Displacement Amplification 2, 3

MLP multilayer perceptron 5, 6, 18

ONT Oxford Nanopore Technologies 4, 9, 10, 15, 16

1059 **PTA** Primary Template-directed Amplification [2](#)

1060

1061 **SA** Supplementary Alignment [16](#), [19](#), [21](#)

1062 **SV** Structural Variation [1–5](#), [9–16](#), [20](#)

1063

1064 **TRA** translocation [9](#), [11](#)

1065

1066 **WGA** Whole Genome Amplification [1–22](#)

1067

1068 **References**

1069

1070 [1] Kalef-Ezra, E. *et al.* Single-cell somatic copy number variants in brain using
1071 different amplification methods and reference genomes. *Communications Biology*
1072 1288 (2024).

1073

1074 [2] Sun, C. *et al.* Mapping recurrent mosaic copy number variation in human neurons.
1075 *Nature Communications* 4220 (2024).

1076

1077 [3] Navin, N. *et al.* Tumour evolution inferred by single-cell sequencing. *Nature* **472**,
1078 90–94 (2011).

1079

1080 [4] Macaulay, I. C. & Voet, T. Single cell genomics: Advances and future perspectives.
1081 *PLOS Genetics* **10**, e1004126 (2014).

1082

1083 [5] Leung, M. L. *et al.* Highly multiplexed targeted dna sequencing from single nuclei.
1084 *Nature Protocols* 214–235 (2016).

1085

1086 [6] Gawad, C., Koh, W. & Quake, S. R. Single-cell genome sequencing: current state
1087 of the science. *Nature Reviews Genetics* 175–188 (2016).

1088

1089 [7] Chen, C. *et al.* Single-cell whole-genome analyses by linear amplification via
1090 transposon insertion (LIANTI). *Science (new York, N.Y.)* **356**, 189–194 (2017).

1091

1092 [8] Zong, C., Lu, S., Chapman, A. R. & Xie, X. S. Genome-wide detection of single-
1093 nucleotide and copy-number variations of a single human cell. *Science* 1622–1626
(2012).

1094

1095 [9] Huang, L., Ma, F., Chapman, A., Lu, S. & Xie, X. S. Single-cell whole-genome
1096 amplification and sequencing: methodology and applications. *Annual Review of*
1097 *Genomics and Human Genetics* 79–102 (2015).

1098

1099 [10] Dean, F. B. *et al.* Comprehensive human genome amplification using multiple
1100 displacement amplification. *Proceedings of the National Academy of Sciences* **99**,
1101 5261–5266 (2002).

1102

1103 [11] Lasken, R. S. & Stockwell, T. B. Mechanism of chimera formation during the
1104 multiple displacement amplification reaction. *BMC Biotechnology* **7**, 19 (2007).

[12]	Pinard, R. <i>et al.</i> Assessment of whole genome amplification-induced bias through high-throughput, massively parallel whole genome sequencing. <i>BMC Genomics</i> 7 , 216 (2006).	1105 1106 1107 1108
[13]	Telenius, H. <i>et al.</i> Degenerate oligonucleotide-primed PCR: General amplification of target DNA by a single degenerate primer. <i>Genomics</i> 13 , 718–725 (1992).	1109 1110 1111
[14]	Gonzalez-Pena, V. <i>et al.</i> Accurate genomic variant detection in single cells with primary template-directed amplification. <i>Proceedings of the National Academy of Sciences</i> 118 , e2024176118 (2021).	1112 1113 1114
[15]	Hård, J. <i>et al.</i> Long-read whole-genome analysis of human single cells. <i>Nature Communications</i> 14 , 5164 (2023).	1115 1116 1117
[16]	Dippenaar, A. <i>et al.</i> Droplet based whole genome amplification for sequencing minute amounts of purified mycobacterium tuberculosis DNA. <i>Scientific Reports</i> 14 , 9931 (2024).	1118 1119 1120 1121
[17]	de Bourcy, C. F. A. <i>et al.</i> A quantitative comparison of single-cell whole genome amplification methods. <i>PLoS ONE</i> e105585 (2014).	1122 1123 1124
[18]	Biezuner, T. <i>et al.</i> Comparison of seven single cell whole genome amplification commercial kits using targeted sequencing. <i>Scientific Reports</i> 17171 (2021).	1125 1126 1127
[19]	Fu, Y. <i>et al.</i> Uniform and accurate single-cell sequencing based on emulsion whole-genome amplification. <i>Proceedings of the National Academy of Sciences</i> 11923–11928 (2015).	1128 1129 1130 1131
[20]	Agyabeng-Dadzie, F. <i>et al.</i> Evaluating the benefits and limits of multiple displacement amplification with whole-genome oxford nanopore sequencing. <i>Molecular Ecology Resources</i> e14094 (2025).	1132 1133 1134
[21]	Dean, F. B., Nelson, J. R., Giesler, T. L. & Lasken, R. S. Rapid amplification of plasmid and phage DNA using Phi29 DNA polymerase and multiply-primed rolling circle amplification. <i>Genome Research</i> 11 , 1095–1099 (2001).	1135 1136 1137 1138
[22]	Lu, N., Qiao, Y., Lu, Z. & Tu, J. Chimera: The spoiler in multiple displacement amplification. <i>Computational and Structural Biotechnology Journal</i> 1688–1696 (2023).	1139 1140 1141 1142
[23]	Lu, N. <i>et al.</i> Exploration of whole genome amplification generated chimeric sequences in long-read sequencing data. <i>Briefings in Bioinformatics</i> 24 , bbad275 (2023).	1143 1144 1145 1146
[24]	Sedlazeck, F. J. <i>et al.</i> Accurate detection of complex structural variations using single-molecule sequencing. <i>Nature Methods</i> 461–468 (2018).	1147 1148 1149 1150

- 1151 [25] Smolka, M. *et al.* Detection of mosaic and population-level structural variants
1152 with sniffles2. *Nature Biotechnology* 1571–1580 (2024).
1153
- 1154 [26] Chen, Y. *et al.* Deciphering the exact breakpoints of structural variations using
1155 long sequencing reads with DeBreak. *Nature Communications* 283 (2023).
1156
- 1157 [27] Heller, D. & Vingron, M. SVIM: Structural variant identification using mapped
1158 long reads. *Bioinformatics* 2907–2915 (2019).
1159
- 1160 [28] Jiang, T. *et al.* Long-read-based human genomic structural variation detection
1161 with cuteSV. *Genome Biology* 189 (2020).
1162
- 1163 [29] Alkan, C., Coe, B. P. & Eichler, E. E. Genome structural variation discovery and
1164 genotyping. *Nature Reviews Genetics* **12**, 363–376 (2011).
1165
- 1166 [30] Kiguchi, Y., Nishijima, S., Kumar, N., Hattori, M. & Suda, W. Long-read
1167 metagenomics of multiple displacement amplified DNA of low-biomass human gut
1168 phageomes by SACRA pre-processing chimeric reads. *DNA Research* **28**, dsab019
1169 (2021).
1170
- 1171 [31] Kosugi, S. *et al.* Comprehensive evaluation of structural variation detection
1172 algorithms for whole genome sequencing. *Genome Biology* **20**, 117 (2019).
1173
- 1174 [32] Mahmoud, M. *et al.* Structural variant calling: The long and the short of it.
1175 *Genome Biology* **20**, 246 (2019).
1176
- 1177 [33] Dalla-Torre, H. *et al.* Nucleotide transformer: building and evaluating robust
1178 foundation models for human genomics. *Nature Methods* 287–297 (2025).
1179
- 1180 [34] Zhou, Z. *et al.* DNABERT-2: Efficient foundation model and benchmark for
1181 multi-species genomes, 1–24 (OpenReview.net, 2024).
1182
- 1183 [35] Nguyen, E. *et al.* HyenaDNA: Long-range genomic sequence modeling at single
1184 nucleotide resolution, Vol. 36, 43177–43201 (Curran Associates, Inc., 2023).
1185
- 1186 [36] Consens, M. E. *et al.* To transformers and beyond: Large language models for
1187 the genome (2023). [arXiv:2311.07621](https://arxiv.org/abs/2311.07621).
1188
- 1189 [37] Routhier, E. & Mozziconacci, J. Genomics enters the deep learning era. *PeerJ*
1190 **10**, e13613 (2022).
1191
- 1192 [38] Li, Y. *et al.* A genomic language model for chimera artifact detection in nanopore
1193 direct rna sequencing. *bioRxiv* (2024). URL <https://www.biorxiv.org/content/early/2024/10/25/2024.10.23.619929>.
1194
- 1195 [39] Poli, M. *et al.* Hyena hierarchy: Towards larger convolutional language models,
1196 Vol. 202, 28043–28078 (PMLR, 2023).

[40] PLC., O. N. Dorado. https://github.com/nanoporetech/dorado (2023).	1197
	1198
[41] Martin, M. Cutadapt removes adapter sequences from high-throughput sequencing reads. <i>Embnnet.journal</i> 17 , 10–12 (2011).	1199
	1200
	1201
[42] Li, H. Minimap2: Pairwise alignment for nucleotide sequences. <i>Bioinformatics</i> 3094–3100 (2018).	1202
	1203
	1204
[43] Danecek, P. <i>et al.</i> Twelve years of SAMtools and BCFtools. <i>GigaScience</i> giab008 (2021).	1205
	1206
	1207
[44] De Coster, W. & Rademakers, R. NanoPack2: Population-scale evaluation of long-read sequencing data. <i>Bioinformatics</i> 39 , btad311 (2023).	1208
	1209
[45] Paszke, A. <i>et al.</i> <i>PyTorch: An imperative style, high-performance deep learning library</i> , Vol. 32, 8024–8035 (Curran Associates, Inc., 2019).	1210
	1211
	1212
[46] Falcon, W. & The PyTorch Lightning team. PyTorch Lightning. GitHub repository (2019). URL https://github.com/Lightning-AI/lightning .	1213
	1214
	1215
[47] Loshchilov, I. & Hutter, F. <i>Decoupled weight decay regularization</i> (2019).	1216
	1217
[48] Yadan, O. Hydra - a framework for elegantly configuring complex applications. GitHub repository (2019). URL https://github.com/facebookresearch/hydra .	1218
	1219
	1220
[49] Chen, X. <i>et al.</i> Manta: Rapid detection of structural variants and indels for germline and cancer sequencing applications. <i>Bioinformatics</i> 1220–1222 (2016).	1221
	1222
	1223
[50] Rausch, T. <i>et al.</i> DELLY: Structural variant discovery by integrated paired-end and split-read analysis. <i>Bioinformatics</i> i333–i339 (2012).	1224
	1225
	1226
[51] Wala, J. A. <i>et al.</i> SvABA: Genome-wide detection of structural variants and indels by local assembly. <i>Genome Research</i> 581–591 (2018).	1227
	1228
	1229
[52] Guo, Q., Li, Y., Wang, T.-Y., Ramakrishnan, A. & Yang, R. OctopusSV and TentacleSV: A one-stop toolkit for multi-sample, cross-platform structural variant comparison and analysis. <i>Bioinformatics</i> btaf599 (2025).	1230
	1231
	1232
[53] English, A. C., Menon, V. K., Gibbs, R. A., Metcalf, G. A. & Sedlazeck, F. J. Truvari: Refined structural variant comparison preserves allelic diversity. <i>Genome Biology</i> 23 , 271 (2022).	1233
	1234
	1235
	1236
[54] Virtanen, P. <i>et al.</i> SciPy 1.0: Fundamental algorithms for scientific computing in python. <i>Nature Methods</i> 261–272 (2020).	1237
	1238
	1239
[55] Hunter, J. D. Matplotlib: A 2d graphics environment. <i>Computing in Science & Engineering</i> 90–95 (2007).	1240
	1241
	1242

1243 [56] Waskom, M. L. seaborn: statistical data visualization. *Journal of Open Source*
1244 *Software* 3021 (2021).

1245
1246
1247
1248
1249
1250
1251
1252
1253
1254
1255
1256
1257
1258
1259
1260
1261
1262
1263
1264
1265
1266
1267
1268
1269
1270
1271
1272
1273
1274
1275
1276
1277
1278
1279
1280
1281
1282
1283
1284
1285
1286
1287
1288